

兔抗 NOC3L 多克隆抗体

中文名称: 兔抗 NOC3L 多克隆抗体

英文名称: Anti-NOC3L rabbit polyclonal antibody

别名: AD24; FAD24; C10orf117

储存: 冷冻 (-20℃) 避光

抗原: NOC3L

宿主: Rabbit

反应种属: Human Mouse

标记物: Unconjugate

相关类别: 一抗

技术规格

Background:	May be required for adipogenesis.
Applications:	WB
Name of antibody:	NOC3L
Immunogen:	Synthesized peptide derived from internal of human NOC3L.
Full name:	NOC3 like DNA replication regulator
Synonyms:	AD24; FAD24; C10orf117
SwissProt:	Q8WTT2
WB Predicted band size:	93 kDa
WB Positive control:	K562 cells lysate
WB Recommended dilution:	500-3000

