

兔抗 RFLNA 多克隆抗体

中文名称：兔抗 RFLNA 多克隆抗体

英文名称：Anti-RFLNA rabbit polyclonal antibody

别名：refilin A; CFM2; FAM101A

相关类别：一抗

储存：冷冻（-20℃）

抗原：RFLNA

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Involved in the regulation of the perinuclear actin network and nuclear shape through interaction with filamins. Plays an essential role in actin cytoskeleton formation in developing cartilaginous cells.
Applications:	ELISA, IHC
Name of antibody:	RFLNA
Immunogen:	Fusion protein of human RFLNA
Full name:	refilin A
Synonyms:	CFM2; FAM101A
SwissProt:	Q6ZTI6

IHC positive control:	Human prostate cancer and Human esophagus cancer
IHC Recommend dilution:	25-100

