

## Alexa Fluor 488 标记山羊抗大鼠 IgG

中文名称: Alexa Fluor 488 标记山羊抗大鼠 IgG

英文名称: Alexa Fluor 488-conjugated Goat anti-rat IgG

相关类别: 荧光标记二抗

储 存: 冷冻 (-20℃)

标 记 物: Alexa Fluor 488

### 技术规格

<b>Background</b>	Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule rat IgG. It also reacts with the light chains of other rat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.
<b>Applications</b>	1/50-200 for most application