

## 兔抗 NOL9 多克隆抗体

中文名称：兔抗 NOL9 多克隆抗体

英文名称： Anti-NOL9 rabbit polyclonal antibody

别名： Grc3; NET6

相关类别： 一抗

抗原： NOL9

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

|                          |  |
|--------------------------|--|
| <b>Background:</b>       | Nucleolar protein 9 is a protein that in humans is encoded by the NOL9 gene. Polynucleotide 5'-kinase involved in rRNA processing. The kinase activity is required for the processing of the 32S precursor into 5.8S and 28S rRNAs, more specifically for the generation of the major 5.8S(S) form. In vitro, has both DNA and RNA 5'-kinase activities. Probably binds RNA. |
| <b>Applications:</b>     | ELISA, WB  |
| <b>Name of antibody:</b> | NOL9   |
| <b>Immunogen:</b>        | Fusion protein of human NOL9   |
| <b>Full name:</b>        | nucleolar protein 9  |

|                                    |            |
|------------------------------------|------------|
| <b>Synonyms :</b>                  | Grc3; NET6 |
| <b>SwissProt:</b>                  | Q5SY16     |
| <b>ELISA Recommended dilution:</b> | 2000-5000  |
| <b>WB Predicted band size:</b>     | 79 kDa     |
| <b>WB Positive control:</b>        | Hela cells |
| <b>WB Recommended dilution:</b>    | 500-2000   |

