

兔抗 NOLC1 多克隆抗体

中文名称: 兔抗 NOLC1 多克隆抗体

英文名称: Anti-NOLC1 rabbit polyclonal antibody

别名: nucleolar and coiled-body phosphoprotein 1; P130; NOPP130; NOPP140; NS5ATP13

储存: 冷冻 (-20°C)

抗原: NOLC1

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Related to nucleologenesis, may play a role in the maintenance of the fundamental structure of the fibrillar center and dense fibrillar component in the nucleolus. It has intrinsic GTPase and ATPase activities. May play an important role in transcription catalyzed by RNA polymerase I.
Applications:	ELISA, IHC
Name of antibody:	NOLC1
Immunogen:	Fusion protein of human NOLC1
Full name:	nucleolar and coiled-body phosphoprotein 1
Synonyms:	P130; NOPP130; NOPP140; NS5ATP13
SwissProt:	Q14978

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-100

