

## 兔抗 RGP1 多克隆抗体

中文名称：兔抗 RGP1 多克隆抗体

英文名称： Anti-RGP1 rabbit polyclonal antibody

别名： RGP1 homolog, RAB6A GEF complex partner 1; KIAA0258

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RGP1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The RIC1-RGP1 complex acts as a guanine nucleotide exchange factor (GEF), which activates RAB6A by exchanging bound GDP for free GTP and may thereby be required for efficient fusion of endosome-derived vesicles with the Golgi compartment. The RIC1-RGP1 complex participates in the recycling of mannose-6-phosphate receptors.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RGP1
<b>Immunogen:</b>	Fusion protein of human RGP1
<b>Full name:</b>	RGP1 homolog, RAB6A GEF complex partner 1

<b>Synonyms:</b>	KIAA0258
<b>SwissProt:</b>	Q92546
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	30-150

