

## 兔抗 NOP58 多克隆抗体

中文名称：兔抗 NOP58 多克隆抗体

英文名称：Anti-NOP58 rabbit polyclonal antibody

别名：NOP58 ribonucleoprotein; NOP5; HSPC120; NOP5/NOP58

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NOP58

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

<b>Background:</b>	The protein encoded by this gene is a core component of box C/D small nucleolar ribonucleoproteins. Some box C/D small nucleolar RNAs (snoRNAs), such as U3, U8, and U14, are dependent upon the encoded protein for their synthesis. This protein is SUMOylated, which is necessary for high affinity binding to snoRNAs.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	NOP58
<b>Immunogen:</b>	Fusion protein of human NOP58
<b>Full name:</b>	NOP58 ribonucleoprotein
<b>Synonyms:</b>	NOP5; HSPC120; NOP5/NOP58

<b>SwissProt:</b>	Q9Y2X3
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human breast cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	59 kDa
<b>WB Positive control:</b>	K562, 231, HeLa, HepG2, A172, 293T cell lysates
<b>WB Recommended dilution:</b>	500-2000



