

兔抗 UBE2F 多克隆抗体

中文名称：兔抗 UBE2F 多克隆抗体

英文名称： Anti-UBE2F rabbit polyclonal antibody

别名： ubiquitin conjugating enzyme E2 F (putative); NCE2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： UBE2F

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Accepts the ubiquitin-like protein NEDD8 from the UBA3-NAE1 E1 complex and catalyzes its covalent attachment to other proteins. The specific interaction with the E3 ubiquitin ligase RBX2, but not RBX1, suggests that the RBX2-UBE2F complex neddylates specific target proteins, such as CUL5.
Applications:	ELISA, IHC
Name of antibody:	UBE2F
Immunogen:	Fusion protein of human UBE2F
Full name:	ubiquitin conjugating enzyme E2 F (putative)
Synonyms:	NCE2

SwissProt:	Q969M7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	50-200

