

Alexa Fluor 555 标记山羊抗小鼠 IgG

中文名称: Alexa Fluor 555 标记山羊抗小鼠 IgG

英文名称: Alexa Fluor 555-conjugated Goat anti-mouse IgG

相关类别: 荧光标记二抗

储 存: 冷冻(-20℃)

标 记 物: VAlexa Fluor 555

技术规格

Background	Immunoglobulin G (IgG) is one of the most abundant protein s in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defense against microorganis ms and the molecules are produced by B lymphocytes as a p art of our adaptive immune response. The IgG molecule has t wo separate functions, to bind to the pathogen that elicited t he response and to recruit other cells and molecules to destr oy the antigen. The variability of the IgG pool is generated b y somatic recombination and the number of specificities in a n individual at a given time point is estimated to be 1011 va riants.
Applications	1/50-200 for most application