

兔抗 NOS2(Ab-151) 多克隆抗体

中文名称: 兔抗 NOS2(Ab-151) 多克隆抗体

英文名称: Anti-NOS2(Ab-151) rabbit polyclonal antibody

别名: NOS; INOS; NOS2A; HEP-NOS

抗原: NOS2(Ab-151)

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

反应种属: Human Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	Produces nitric oxide (NO) which is a messenger molecule with diverse functions throughout the body. In macrophages, NO mediates tumoricidal and bactericidal actions. Also has nitrosylase activity and mediates cysteine S-nitrosylation of cytoplasmic target proteins such COX2.
Applications:	WB
Name of antibody:	NOS2(Ab-151)
Immunogen:	Synthesized non-phosphopeptide derived from human iNOS around the phosphorylation site of tyrosine 151 (Q-Y-Y(p)-G-S).
Full name:	nitric oxide synthase 2, inducible

Synonyms :	NOS; INOS; NOS2A; HEP-NOS
SwissProt:	P35228
WB Predicted band size:	131 kDa
WB Positive control:	NIH/3T3 cells lysate
WB Recommended dilution:	500-3000

