

## 兔抗 IGFBP5 多克隆抗体

中文名称:	兔抗 IGFBP5 多克隆抗体
英文名称:	Anti-IGFBP5 rabbit polyclonal antibody
别名:	IBP5
抗原:	IGFBP5
储存:	冷冻 (-20℃)
宿主:	Rabbit
相关类别:	一抗
反应种属:	Human, Mouse, Rat
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

### 技术规格

<b>Background:</b>	Insulin-like growth factor-binding protein 5 is a protein that in humans is encoded by the IGFBP5 gene. IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	IGFBP5
<b>Immunogen:</b>	Fusion protein of human IGFBP5
<b>Full name:</b>	insulin-like growth factor binding protein 5
<b>Synonyms :</b>	IBP5

<b>SwissProt:</b>	P24593
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human liver cancer and Human thyroid cancer
<b>IHC Recommend dilution :</b>	25-100

