

兔抗 RGS10 多克隆抗体

中文名称：兔抗 RGS10 多克隆抗体

英文名称： Anti-RGS10 rabbit polyclonal antibody

储 存： 冷冻（-20℃）

抗 原： RGS10

宿 主： Rabbit

反应种属： Human, Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Regulator of G protein signaling (RGS) family members are regulatory molecules that act as GTPase activating proteins (GAPs) for G alpha subunits of heterotrimeric G proteins. RGS proteins are able to deactivate G protein subunits of the Gi alpha, Go alpha and Gq alpha subtypes. They drive G proteins into their inactive GDP-bound forms. Regulator of G protein signaling 10 belongs to this family. All RGS proteins share a conserved 120-amino acid sequence termed the RGS domain.
Applications:	ELISA, IHC
Name of antibody:	RGS10
Immunogen:	Fusion protein of human RGS10

Full name:	regulator of G-protein signaling 10
SwissProt:	O43665
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	10-50

