

兔抗 NOS2 多克隆抗体

中文名称：兔抗 NOS2 多克隆抗体

英文名称： Anti-NOS2 rabbit polyclonal antibody

别名： NOS; INOS; NOS2A; HEP-NOS

抗原： NOS2

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Nitric oxide is a reactive free radical which acts as a biologic mediator in several processes, including neurotransmission and antimicrobial and antitumoral activities. This gene encodes a nitric oxide synthase which is expressed in liver and is inducible by a combination of lipopolysaccharide and certain cytokines. Three related pseudogenes are located within the Smith-Magenis syndrome region on chromosome 17.
Applications:	ELISA, IHC
Name of antibody:	NOS2
Immunogen:	Synthetic peptide of human NOS2
Full name:	nitric oxide synthase 2, inducible

Synonyms :	NOS; INOS; NOS2A; HEP-NOS
SwissProt:	P35228
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and Human lung cancer
IHC Recommend dilution:	25-100

