

兔抗 NOS3 多克隆抗体

中文名称：兔抗 NOS3 多克隆抗体

英文名称：Anti-NOS3 rabbit polyclonal antibody

级 别：eNOS; ECNOS

抗 原：NOS3

储 存：冷冻（-20℃）

宿 主：Rabbit

反应种属：Human

相关类别：一抗

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Nitric oxide is a reactive free radical which acts as a biologic mediator in several processes, including neurotransmission and antimicrobial and antitumoral activities. Nitric oxide is synthesized from L-arginine by nitric oxide synthases. Variations in this gene are associated with susceptibility to coronary spasm. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	NOS3
Immunogen:	Synthetic peptide of human NOS3
Full name:	nitric oxide synthase 3

Synonyms:	eNOS; ECNOS
SwissProt:	P29474
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	10-50



