

## 兔抗 NOSIP 多克隆抗体

中文名称：兔抗 NOSIP 多克隆抗体

英文名称： Anti-NOSIP rabbit polyclonal antibody

别名： CGI-25

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NOSIP

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

|                          |  |
|--------------------------|--|
| <b>Background:</b>       | The protein encoded by this gene may modulate the activity and localization of nitric oxide synthase (endothelial and neuronal) and thus nitric oxide production. Alternative splicing results in multiple transcript variants that encode the same protein. |
| <b>Applications:</b>     | ELISA, IHC   |
| <b>Name of antibody:</b> | NOSIP  |
| <b>Immunogen:</b>        | Full length fusion protein   |
| <b>Full name:</b>        | nitric oxide synthase interacting protein  |
| <b>Synonyms:</b>         | CGI-25   |
| <b>SwissProt:</b>        | Q9Y314   |

|                                    |   |
|------------------------------------|---|
| <b>ELISA Recommended dilution:</b> | 5000-10000                              |
| <b>IHC positive control:</b>       | Human liver cancer;Human gastric cancer |
| <b>IHC Recommend dilution:</b>     | 25-100                                  |



