

## KIRREL2 抗原(重组蛋白)

中文名称: KIRREL2 抗原(重组蛋白)

英文名称: KIRREL2 Antigen (Recombinant Protein)

别 名: kirre like nephrin family adhesion molecule 2; NLG1; NEPH3; FILTRIN

储 存: 冷冻(-20℃)

相关类别: 抗原

## 概述

Fusion protein corresponding to C terminal 170 amino acids of human KIRREL2

## 技术规格

Full name:	kirre like nephrin family adhesion molecule 2
Synonyms:	NLG1; NEPH3; FILTRIN
Swissprot:	Q6UWL6
Gene Accession:	BC007312
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a type I transmembrane protein and member of the immun oglobulin superfamily of cell adhesion molecules. The encoded protein localizes to adherens junctions in pancreatic beta cells and regulates insulin secretion. A utoantibodies against the encoded protein have been detected in serum from patients with type 1 diabetes. This gene may also play a role in glomerular de velopment and decreased expression of this gene has been observed in huma n glomerular diseases. This gene and the related opposite-strand gene nephrin



(GeneID: 527362) are regulated by a bidirectional promoter.