

Cy3 标记小鼠抗大鼠 IgG

- 中文名称: Cy3 标记小鼠抗大鼠 IgG
- 英文名称: Cy3-conjugated Mouse anti-rat IgG
- 相关类别: 荧光标记二抗
- 储 存: 冷冻 (-20℃)
- 标记物: Cy3

技术规格

| Background: | Immunoglobulin G (IgG) is one of the most abundant proteins in ser um with normal levels between 8-17 mg/ml in adult blood. IgG is im portant for our defense against microorganisms and the molecules ar e produced by B lymphocytes as a part of our adaptive immune resp onse. The IgG molecule has two separate functions; to bind to the p athogen that elicited the response and to recruit other cells and mol ecules to destroy the antigen. The variability of the IgG pool is gener ated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants. |
|---------------|--|
| Applications: | 1/100-1000 for Immunofluorescence |