

KLHDC2 抗原（重组蛋白）

中文名称：KLHDC2 抗原（重组蛋白）

英文名称：KLHDC2 Antigen (Recombinant Protein)

别名：kelch domain containing 2; LCP; HCLP1; HCLP-1

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 207-406 amino acids of human KLHDC2

技术规格

Full name:	kelch domain containing 2
Synonyms:	LCP; HCLP1; HCLP-1
Swissprot:	Q9Y2U9
Gene Accession:	BC002335
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Represses CREB3-mediated transcription by interfering with CREB3-DNA binding.