

## KLHL1 抗原（重组蛋白）

中文名称： KLHL1 抗原（重组蛋白）

英文名称： KLHL1 Antigen (Recombinant Protein)

别名： kelch-like family member 1; MRP2

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to C terminal 250 amino acids of human KLHL1

### 技术规格

<b>Full name:</b>	kelch-like family member 1
<b>Synonyms:</b>	MRP2
<b>Swissprot:</b>	Q9NR64
<b>Gene Accession:</b>	BC022460
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The KLHL1 protein belongs to a family of actin-organizing proteins related to Drosophila Kelch (Nemes et al., 2000 [PubMed 10888605]).