

兔抗 IGSF3 多克隆抗体

中文名称：兔抗 IGSF3 多克隆抗体

英文名称：Anti-IGSF3 rabbit polyclonal antibody

别名：immunoglobulin superfamily member 3 V8; LCDD; EWI-3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：IGSF3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is an immunoglobulin-like membrane protein containing several V-type Ig-like domains. A mutation in this gene has been associated with bilateral nasolacrimal duct obstruction (LCDD).
Applications:	ELISA, IHC
Name of antibody:	IGSF3
Immunogen:	Synthetic peptide of human IGSF3
Full name:	immunoglobulin superfamily member 3

Synonyms:	V8; LCDD; EWI-3
SwissProt:	O75054
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer; Human liver cancer
IHC Recommend dilution:	20-100



