

兔抗 UBE2R2 多克隆抗体

中文名称：兔抗 UBE2R2 多克隆抗体

英文名称：Anti-UBE2R2 rabbit polyclonal antibody

别名：ubiquitin conjugating enzyme E2 R2; UBC3B; CDC34B; E2-CDC34B

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：UBE2R2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Protein kinase CK2 is a ubiquitous and pleiotropic Ser/Thr protein kinase involved in cell growth and transformation. This gene encodes a protein similar to the E2 ubiquitin conjugating enzyme UBC3/CDC34. Studies suggest that CK2-dependent phosphorylation of this ubiquitin-conjugating enzyme functions by regulating beta-TrCP substrate recognition and induces its interaction with beta-TrCP, enhancing beta-catenin degradation.
Applications:	ELISA, IHC
Name of antibody:	UBE2R2

Immunogen:	Fusion protein of human UBE2R2
Full name:	ubiquitin conjugating enzyme E2 R2
Synonyms:	UBC3B; CDC34B; E2-CDC34B
SwissProt:	Q712K3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer and Human tonsil
IHC Recommend dilution:	150-300



