

兔抗 UBE2S 多克隆抗体

中文名称：兔抗 UBE2S 多克隆抗体

英文名称： Anti-UBE2S rabbit polyclonal antibody

别名： EPF5; E2EPF; E2-EPF

相关类别： 一抗

抗原： UBE2S

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the ubiquitin-conjugating enzyme family. The encoded protein is able to form a thiol ester linkage with ubiquitin in a ubiquitin activating enzyme-dependent manner, a characteristic property of ubiquitin carrier proteins.
Applications:	ELISA, WB, IHC
Name of antibody:	UBE2S
Immunogen:	Fusion protein of human UBE2S
Full name:	ubiquitin-conjugating enzyme E2S
Synonyms :	EPF5; E2EPF; E2-EPF
SwissProt:	Q16763

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human brain
IHC Recommend dilution:	50-200
WB Predicted band size:	24 kDa
WB Positive control:	Human testis tissue
WB Recommended dilution:	200-1000



