

兔抗 NOTCH3 多克隆抗体

中文名称：兔抗 NOTCH3 多克隆抗体

英文名称：Anti-NOTCH3 rabbit polyclonal antibody

别名：CASIL, CADASIL

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NOTCH3

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes the third discovered human homologue of the *Drosophila melanogaster* type I membrane protein notch. In *Drosophila*, notch interaction with its cell-bound ligands (delta, serrate) establishes an intercellular signalling pathway that plays a key role in neural development. Homologues of the notch-ligands have also been identified in human, but precise interactions between these ligands and the human notch homologues remains to be determined. Mutations in NOTCH3 have been identified as the underlying cause of cerebral autosomal dominant

	ant arteriopathy with subcortical infarcts and leukoencephalopathy (CADASIL).
Applications:	ELISA, WB, IHC
Name of antibody:	NOTCH3
Immunogen:	Synthetic peptide of human NOTCH3
Full name:	Notch 3
Synonyms :	CASIL, CADASIL
SwissProt:	Q9UM47
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human colon cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	244 kDa
WB Positive control:	Human colon cancer tissue
WB Recommended dilution:	200-1000



