

## 兔抗 RHOBTB1 多克隆抗体

中文名称：兔抗 RHOBTB1 多克隆抗体

英文名称： Anti-RHOBTB1 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： RHOBTB1

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene belongs to the Rho family of the small GTPase superfamily. It contains a GTPase domain, a proline-rich region, a tandem of 2 BTB (broad complex, tramtrack, and bric-a-brac) domains, and a conserved C-terminal region. The protein plays a role in small GTPase-mediated signal transduction and the organization of the actin filament system.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	RHOBTB1
<b>Immunogen:</b>	Fusion protein of human RHOBTB1
<b>Full name:</b>	Rho-related BTB domain containing 1
<b>SwissProt:</b>	O94844
<b>ELISA Recommended dilution:</b>	2000-5000

<b>WB Predicted band size:</b>	79 kDa
<b>WB Positive control:</b>	A431 cells
<b>WB Recommended dilution:</b>	500-2000

