

兔抗 RHOBTB2 多克隆抗体

中文名称：兔抗 RHOBTB2 多克隆抗体

英文名称：Anti-RHOBTB2 rabbit polyclonal antibody

别名：DBC2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RHOBTB2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a small Rho GTPase and a candidate tumor suppressor. The encoded protein interacts with the cullin-3 protein, a ubiquitin E3 ligase necessary for mitotic cell division. This protein inhibits the growth and spread of some types of breast cancer. Three transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	RHOBTB2
Immunogen:	Synthetic peptide of human RHOBTB2
Full name:	Rho-related BTB domain containing 2

Synonyms :	DBC2
SwissProt:	Q9BYZ6
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and Human gastric cancer
IHC Recommend dilution:	25-100

