

## 兔抗 UBE2U 多克隆抗体

中文名称：兔抗 UBE2U 多克隆抗体

英文名称：Anti-UBE2U rabbit polyclonal antibody

抗 原： UBE2U

储 存： 冷冻（-20℃）

宿 主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	UBE2U (ubiquitin-conjugating enzyme E2U) is a 321 amino acid E2 ubiquitin conjugating enzyme that catalyzes the covalent attachment of ubiquitin to other proteins. Existing as two alternatively spliced isoforms, UBE2U is encoded by a gene located on human chromosome 1, which spans 260 million base pairs, contains over 3,000 genes and comprises nearly 8% of the human genome.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UBE2U
<b>Immunogen:</b>	Fusion protein of human UBE2U
<b>Full name:</b>	ubiquitin-conjugating enzyme E2U (putative)
<b>SwissProt:</b>	Q5VVX9

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human breast cancer
<b>IHC Recommend dilution:</b>	25-100

