

兔抗 NOX5 多克隆抗体

中文名称：兔抗 NOX5 多克隆抗体

英文名称：Anti-NOX5 rabbit ployclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：NOX5

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

WB Recommended dilution:	200-1000
WB Positive control:	Raji, TM4 and PC3 cells
WB Predicted band size:	86 kDa
ELISA Recommended dilution:	1000-2000
SwissProt:	Q96PH1
Full name:	NADPH oxidase, EF-hand calcium binding domain 5
Immunogen:	Fusion protein of human NOX5
Name of antibody:	NOX5
Applications:	ELISA, WB
Background:	This gene is predominantly expressed in the testis and lymphocyte-rich areas of spleen and lymph nodes. It encodes a calcium-dependen NADPH oxi

dase that generates superoxide, and functions as a calcium-dependent proton channel that may regulate redox-dependent processes in lymphocytes and spermatozoa. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

