

FITC 标记驴抗小鼠 IgG

中文名称： FITC 标记驴抗小鼠 IgG

英文名称： FITC-conjugated Donkey Anti-Mouse IgG

相关类别： 荧光标记二抗

储 存： 冷冻 (-20℃)

标 记物： FITC 储存冷冻 (-20℃)

标 记物： FITC

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule mouse IgG. It also reacts with the light chains of other mouse immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with bovine, chicken, goat, guinea pig, Syrian hamster, horse, human, rabbit, rat, and sheep serum proteins, but it may cross-react with immunoglobulins from other species.
Applications	1/50-200 for most application