

兔抗 NPAP1 多克隆抗体

中文名称：兔抗 NPAP1 多克隆抗体

英文名称：Anti-NPAP1 rabbit polyclonal antibody

别名：C15orf2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NPAP1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is located in the Prader-Willi syndrome region on chromosome 15. This gene is biallelically expressed in adult testis and brain but is paternally imprinted in fetal brain. Defects in this gene may be associated with Prader-Willi syndrome.
Applications:	ELISA, IHC
Name of antibody:	NPAP1
Immunogen:	Synthetic peptide of human NPAP1
Full name:	nuclear pore associated protein 1
Synonyms :	C15orf2
SwissProt:	Q9NZP6

ELISA Recommended dilution:	1000-5000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-200

