

兔抗 RHOH 多克隆抗体

中文名称：兔抗 RHOH 多克隆抗体

英文名称：Anti-RHOH rabbit polyclonal antibody

别名：ras homolog family member H; TTF; ARHH

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RHOH

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a member of the Ras superfamily of guanosine triphosphate (GTP)-metabolizing enzymes. The encoded protein is expressed in hematopoietic cells, where it functions as a negative regulator of cell growth and survival. This gene may be hypermutated or misexpressed in leukemias and lymphomas. Chromosomal translocations in non-Hodgkin's lymphoma occur between this locus and B-cell CLL/lymphoma 6 (BCL6) on chromosome 3, leading to the production of fusion transcripts. Alternative splicing in the 5' untranslated region results in multi

	ple transcript variants that encode the same protein.
Applications:	ELISA, IHC
Name of antibody:	RHOH
Immunogen:	Full length fusion protein
Full name:	ras homolog family member H
Synonyms:	TTF; ARHH
SwissProt:	Q15669
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	40-200

