

LCMT1 抗原（重组蛋白）

中文名称：LCMT1 抗原（重组蛋白）

英文名称：LCMT1 Antigen (Recombinant Protein)

别名：leucine carboxyl methyltransferase 1; LCMT; PPMT1; CGI-68

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to C terminal 209 amino acids of human LCMT1

技术规格

Full name:	leucine carboxyl methyltransferase 1
Synonyms:	LCMT; PPMT1; CGI-68
Swissprot:	Q9UIC8
Gene Accession:	BC014217
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	LCMT1 catalyzes the methylation of the carboxyl group of the C-terminal leucine residue (leu309) of the catalytic subunit of protein phosphatase-2A.