

## 兔抗 UBE4A 多克隆抗体

- 中文名称：兔抗 UBE4A 多克隆抗体
- 英文名称：Anti-UBE4A rabbit polyclonal antibody
- 别名：E4; UFD2; UBOX2
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：UBE4A
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the U-box ubiquitin ligase family. The encoded protein is involved in mult iubiquitin chain assembly and plays a critical role in c hromosome condensation and separation through the polyubiquitination of securin. Autoantibodies against t he encoded protein may be markers for scleroderma and Crohn's disease. A pseudogene of this gene is lo cated on the long arm of chromosome 3. Alternativel y spliced transcript variants encoding multiple isoform s have been observed for this gene.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	UBE4A
<b>Immunogen:</b>	Fusion protein of human UBE4A
<b>Full name:</b>	ubiquitination factor E4A
<b>Synonyms:</b>	E4; UFD2; UBOX2
<b>SwissProt:</b>	Q14139
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and human breast cancer
<b>IHC Recommend dilution:</b>	25-100



