

## LILRA2 抗原（重组蛋白）

中文名称：LILRA2 抗原（重组蛋白）

英文名称：LILRA2 Antigen (Recombinant Protein)

别名：leukocyte immunoglobulin like receptor A2; ILT1; LIR7; CD85H; LIR-7

储存：冷冻（-20℃）

相关类别：抗原

### 概述

Fusion protein corresponding to a region derived from 24-223 amino acids of human LILRA2

### 技术规格

<b>Full name:</b>	leukocyte immunoglobulin like receptor A2
<b>Synonyms:</b>	ILT1; LIR7; CD85H; LIR-7
<b>Swissprot:</b>	Q8N149
<b>Gene Accession:</b>	BC017412
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of a family of immunoreceptors that are expressed predominantly on monocytes and B cells, and at lower levels on dendritic cells and natural killer cells. The encoded protein is an activating receptor that inhibits dendritic cell differentiation and antigen presentation and suppresses innate immune response. Alternatively spliced transcript variants encoding different isoforms have been found. This gene is located in a cluster of related

genes on chromosome 19 and there is a pseudogene for this gene on chromosome 3.