

兔抗 IKBKG(Ab-85) 多克隆抗体

中文名称:	兔抗 IKBKG(Ab-85) 多克隆抗体
英文名称:	Anti-ICKBKG(Ab-85) rabbit polyclonal antibody
别名:	IP; IP1; IP2; FIP3; IPD2; NEMO; FIP-3; Fip3p; IMD33; AMCBX1; ZC2HC9; IKK-gamma
相关类别:	一抗
储存:	冷冻 (-20℃) 避光
宿主:	Rabbit
抗原:	IKBKG(Ab-85)
反应种属:	Human
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	This gene encodes the regulatory subunit of the inhibitor of kappaB kinase (IKK) complex, which activates NF-kappaB resulting in activation of genes involved in inflammation, immunity, cell survival, and other pathways. Mutations in this gene result in incontinentia pigmenti, hypohidrotic ectodermal dysplasia, and several other types of immunodeficiencies. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	WB, IF
Name of antibody:	IKBKG(Ab-85)
Immunogen:	Synthesized non-phosphopeptide derived from human IKK- γ

	around the phosphorylation site of serine 85 (Q-A-S(p)-Q-R).
Full name:	inhibitor of kappa light polypeptide gene enhancer in B-cells , kinase gamma
Synonyms :	IP; IP1; IP2; FIP3; IPD2; NEMO; FIP-3; Fip3p; IMD33; AMCBX1; ZC2HC9; IKK-gamma
SwissProt:	Q9Y6K9
WB Predicted band size:	48 kDa
WB Positive control:	HepG2 cells treated with Anisomycin (0.5uM, 5hours) lysate
WB Recommended dilution :	500-3000
IF positive control:	Hela cells
IF Recommend dilution:	100-500



