

FITC 标记兔抗人 IgG

中文名称: FITC 标记兔抗人 IgG

英文名称: FITC-conjugated Rabbit anti-human IgG

相关类别: 荧光标记二抗

储 存: 冷冻 (-20℃)

标 记 物: FITC

技术规格

Background	Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule human IgG. It also reacts with the light chains of other human immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.
Applications	1/50-200 for most application

mlbio 梅联生物
Good elisakit producers