

兔抗 UBFD1 多克隆抗体

中文名称：兔抗 UBFD1 多克隆抗体

英文名称：Anti-UBFD1 rabbit polyclonal antibody

别名：UBPH

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：UBFD1

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	May play a role as NF-kappa-B regulator.
Applications:	WB
Name of antibody:	UBFD1
Immunogen:	Synthesized peptide derived from internal of human UBFD1.
Full name:	ubiquitin family domain containing 1
Synonyms :	UBPH
SwissProt:	O14562
WB Predicted band size:	33 kDa
WB Positive control:	K562 cells lysate
WB Recommended dilution:	500-3000

