

## 兔抗 NPDC1 多克隆抗体

中文名称: 兔抗 NPDC1 多克隆抗体

英文名称: Anti-NPDC1 rabbit polyclonal antibody

别名: CAB; CAB-; CAB1; CAB-1; NPDC-1

储存: 冷冻 (-20℃) 避光

抗原: NPDC1

宿主: Rabbit

反应种属: Human Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

<b>Background:</b>	Suppresses oncogenic transformation in neural and non-neural cells and down-regulates neural cell proliferation. Might be involved in transcriptional regulation By similarity.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	NPDC1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human NPDC 1.
<b>Full name:</b>	neural proliferation, differentiation and control, 1
<b>Synonyms :</b>	CAB; CAB-; CAB1; CAB-1; NPDC-1
<b>SwissProt:</b>	Q9NQX5
<b>WB Predicted band size:</b>	35 kDa

<b>WB Positive control:</b>	LOVO cells and HepG2 cells lysates
<b>WB Recommended dilution:</b>	500-3000

