

## 兔抗 RIMS4 多克隆抗体

中文名称：兔抗 RIMS4 多克隆抗体

英文名称：Anti-RIMS4 rabbit polyclonal antibody

别名：RIM4; RIM 4; RIM-4; C20orf190; RIM4gamma; RIM4-gamma

抗原：RIMS4

储存：冷冻（-20℃）避光

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Regulates synaptic membrane exocytosis By similarity.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	RIMS4
<b>Immunogen:</b>	Synthesized peptide derived from internal of human RIMS4.
<b>Full name:</b>	regulating synaptic membrane exocytosis 4
<b>Synonyms :</b>	RIM4; RIM 4; RIM-4; C20orf190; RIM4gamma; RIM4-gamma
<b>SwissProt:</b>	Q9H426
<b>IHC positive control:</b>	Human breast carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	29 kDa
<b>WB Positive control:</b>	Jurkat cells and Hela cells lysates

WB Recommended dilution: 500-3000



