

兔抗 NPHP4 多克隆抗体

中文名称: 兔抗 NPHP4 多克隆抗体

英文名称: Anti-NPHP4 rabbit polyclonal antibody

别名: POC10; SLSN4

宿主: Rabbit

储存: 冷冻 (-20℃) 避光

抗原: NPHP4

反应种属: Human

标记物: Unconjugate

相关类别: 一抗

克隆类型: Unconjugate

技术规格

Background:	This gene encodes a protein involved in renal tubular development and function. This protein interacts with nephrocystin, and belongs to a multifunctional complex that is localized to actin- and microtubule-based structures. Mutations in this gene are associated with nephronophthisis type 4, a renal disease, and with Senior-Loken syndrome type 4, a combination of nephronophthisis and retinitis pigmentosa. Alternative splicing results in multiple transcript variants.
Applications:	IHC
Name of antibody:	NPHP4

Immunogen:	Synthesized peptide derived from internal of human NPH P4.
Full name:	nephronophthisis 4
Synonyms :	POC10; SLSN4
SwissProt:	O75161
IHC positive control:	Human colon carcinoma tissue
IHC Recommend dilution:	50-100

