

## 兔抗 UBIAD1 多克隆抗体

中文名称：兔抗 UBIAD1 多克隆抗体

英文名称：Anti-UBIAD1 rabbit polyclonal antibody

别名：UbiA prenyltransferase domain containing 1; SCCD; TERE1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：UBIAD1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein thought to be involved in cholesterol and phospholipid metabolism. Mutations in this gene are associated with Schnyder crystalline corneal dystrophy.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UBIAD1
<b>Immunogen:</b>	Synthetic peptide of human UBIAD1
<b>Full name:</b>	UbiA prenyltransferase domain containing 1
<b>Synonyms:</b>	SCCD; TERE1
<b>SwissProt:</b>	Q9Y5Z9
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human liver cancer and Human ovarian cancer
<b>IHC Recommend dilution:</b>	30-150



