

兔抗 RIN3 多克隆抗体

中文名称：兔抗 RIN3 多克隆抗体

英文名称：Anti-RIN3 rabbit polyclonal antibody

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：RIN3

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Ras effector protein that functions as a guanine nucleotide exchange (GEF) for RAB5B and RAB31, by exchanging bound GDP for free GTP. Required for normal RAB31 function.
Applications:	WB
Name of antibody:	RIN3
Immunogen:	Synthesized peptide derived from internal of human RIN3.
Full name:	Ras and Rab interactor 3
Synonyms :	0
SwissProt:	Q8TB24
WB Predicted band size:	108 kDa
WB Positive control:	Huvec cells and K562 cells lysates

WB Recommended dilution: 500-3000

