

## 兔抗 NPHS1 多克隆抗体

中文名称：兔抗 NPHS1 多克隆抗体

英文名称：Anti-NPHS1 rabbit polyclonal antibody

别名：CNF; NPHN; nephrin

相关类别：一抗

抗原：NPHS1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

#### Background:

This gene encodes a member of the immunoglobulin family of cell adhesion molecules that functions in the glomerular filtration barrier in the kidney. The gene is primarily expressed in renal tissues, and the protein is a type-1 transmembrane protein found at the slit diaphragm of glomerular podocytes. The slit diaphragm is thought to function as an ultrafilter to exclude albumin and other plasma macromolecules in the formation of urine. Mutations in this gene result in Finnish-type congenital nephrosis 1, characterized by severe proteinuria and loss of the slit diaphragm and foot processes.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NPHS1
<b>Immunogen:</b>	Synthetic peptide of human NPHS1
<b>Full name:</b>	nephrosis 1, congenital, Finnish type (nephrin)
<b>Synonyms :</b>	CNF; NPHN; nephrin
<b>SwissProt:</b>	O60500
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	25-100

