

兔抗 RIOK1 多克隆抗体

中文名称：兔抗 RIOK1 多克隆抗体

英文名称： Anti-RIOK1 rabbit polyclonal antibody

别名： RIO kinase 1; RIO1; AD034; RRP10; bA288G3.1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RIOK1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene competes with pI Cln for inclusion in the protein arginine methyltransferase 5 complex. This complex targets substrates for dimethylation. The encoded protein is essential for the last steps in the maturation of 40S subunits.
Applications:	ELISA, IHC
Name of antibody:	RIOK1
Immunogen:	Fusion protein of human RIOK1
Full name:	RIO kinase 1
Synonyms:	RIO1; AD034; RRP10; bA288G3.1
SwissProt:	Q9BRS2

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	70-350

