

## 兔抗 RIOK3 多克隆抗体

中文名称：兔抗 RIOK3 多克隆抗体

英文名称： Anti-RIOK3 rabbit polyclonal antibody

别名： SUDD

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： RIOK3

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene was identified by the similarity of its product to the <i>Aspergillus nidulans</i> SUDD protein, an extragenic suppressor of the heat-sensitive bimD6 mutation that fails to attach properly to the spindle microtubules at a restrictive temperature. The specific function of this gene has not yet been determined.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	RIOK3
<b>Immunogen:</b>	Synthesized peptide derived from internal of human RIOK3.
<b>Full name:</b>	RIO kinase 3

<b>Synonyms :</b>	SUDD
<b>SwissProt:</b>	O14730
<b>WB Predicted band size:</b>	59 kDa
<b>WB Positive control:</b>	HepG2 cells lysate
<b>WB Recommended dilution:</b>	500-3000

