

## 兔抗 UBL4B 多克隆抗体

中文名称：兔抗 UBL4B 多克隆抗体

英文名称： Anti-UBL4B rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： UBL4B

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

Ubiquitin is a 77 amino acid protein that targets proteins for degradation by the 26S proteasome. Ubiquitin-like proteins are not directly involved in protein degradation, but appear to have many mechanistic similarities with the ubiquitin pathway. UBL4B (Ubiquitin-like protein 4B) is a 174 amino acid cytoplasmic protein that contains a ubiquitin-like domain. UBL4B likely arose by retroposition during mammalian evolution from UBL4A, an X-linked intron-bearing housekeeping gene. While UBL4A is highly conserved, UBL4B has undergone rapid evolution and may have evolved new functions. Expression of UBL4B is restricted to post-meio

	tic germ cells in testis and ovarian tissue, where it likely functions in post-translational protein modification.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UBL4B
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	ubiquitin-like 4B
<b>SwissProt:</b>	Q8N7F7
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human esophagus cancer and human liver cancer
<b>IHC Recommend dilution:</b>	25-100

