

## 兔抗 NPM1(Phospho-Thr199) 多克隆抗体

中文名称：兔抗 NPM1(Phospho-Thr199) 多克隆抗体

英文名称： Anti-NPM1(Phospho-Thr199) rabbit polyclonal antibody

别名： B23; NPM

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： NPM1(Phospho-Thr199)

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

<b>Background:</b>	This gene encodes a phosphoprotein which moves between the nucleus and the cytoplasm. The gene product is thought to be involved in several processes including regulation of the ARF/p53 pathway. A number of genes are fusion partners have been characterized, in particular the anaplastic lymphoma kinase gene on chromosome 2. Mutations in this gene are associated with acute myeloid leukemia. More than a dozen pseudogenes of this gene have been identified. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	WB

<b>Name of antibody:</b>	NPM1(Phospho-Thr199)
<b>Immunogen:</b>	Peptide sequence around phosphorylation site of Threonine 199 (R-D-T(p)-P-A) derived from Human NPM.
<b>Full name:</b>	nucleophosmin (nucleolar phosphoprotein B23, numatrin)
<b>Synonyms :</b>	B23; NPM
<b>SwissProt:</b>	P06748
<b>WB Predicted band size:</b>	33 kDa
<b>WB Positive control:</b>	Mouse heart tissue lysates
<b>WB Recommended dilution:</b>	500-1000

