

## 兔抗 NPPA 多克隆抗体

中文名称：兔抗 NPPA 多克隆抗体

英文名称： Anti-NPPA rabbit polyclonal antibody

别名： ANF; ANP; PND; ATFB6; CDD-ANF

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： NPPA

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene belongs to the natriuretic peptide family. Natriuretic peptides are implicated in the control of extracellular fluid volume and electrolyte homeostasis. This protein is synthesized as a large precursor (containing a signal peptide), which is processed to release a peptide from the N-terminus with similarity to vasoactive peptide, cardiodilatin, and another peptide from the C-terminus with natriuretic-diuretic activity. Mutations in this gene have been associated with atrial fibrillation familial type 6.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	NPPA
<b>Immunogen:</b>	Fusion protein of human NPPA
<b>Full name:</b>	natriuretic peptide A
<b>Synonyms :</b>	ANF; ANP; PND; ATFB6; CDD-ANF
<b>SwissProt:</b>	P01160
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human brain
<b>IHC Recommend dilution:</b>	25-100

