

## 兔抗 UBL7 多克隆抗体

- 中文名称：兔抗 UBL7 多克隆抗体
- 英文名称：Anti-UBL7 rabbit polyclonal antibody
- 别名：ubiquitin like 7; TCBA1; BMSC-UbP
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：UBL7
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

**Background:**

Ubiquitin is an important cellular signal that targets proteins for degradation or regulates their functions. UBL7 (ORF name SB132) encodes ubiquitin-like protein 7 also named bone marrow stromal cell ubiquitin-like protein (BMSC-UbP) capable of binding ubiquitin. BMSC-UbP protein is derived from bone marrow stromal cells containing a ubiquitin-associated (UBA) domain at the C terminus that has been implicated in linking cellular processes and the ubiquitin system. BMSC-UbP mRNA is widely expressed in human multiple tissues and various tumor cell lines and is decreased in BMS

	C stimulated with phorbol myristate acetate (PMA) suggesting a role in regulation of BMSC function.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	UBL7
<b>Immunogen:</b>	Fusion protein of human UBL7
<b>Full name:</b>	ubiquitin like 7
<b>Synonyms:</b>	TCBA1; BMSC-UbP
<b>SwissProt:</b>	Q96S82
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	100-200
<b>WB Predicted band size:</b>	41 kDa
<b>WB Positive control:</b>	Mouse testis tissue, Hela cell lysates
<b>WB Recommended dilution:</b>	500-2000



