

兔抗 RIT1 多克隆抗体

中文名称：兔抗 RIT1 多克隆抗体

英文名称： Anti-RIT1 rabbit polyclonal antibody

别名： NS8; RIT; RIBB; ROC1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： RIT1

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of a subfamily of Ras-related GTPases. The encoded protein is involved in regulating p38 MAPK-dependent signaling cascades related to cellular stress. This protein also cooperates with nerve growth factor to promote neuronal development and regeneration. Alternate splicing results in multiple transcript variants.
Applications:	WB
Name of antibody:	RIT1
Immunogen:	Synthesized peptide derived from human RIT1.
Full name:	Ras-like without CAAX 1

Synonyms :	NS8; RIT; RIBB; ROC1
SwissProt:	Q92963
WB Predicted band size:	25 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

