

## 兔抗 IKZF4 多克隆抗体

中文名称: 兔抗 IKZF4 多克隆抗体

英文名称: Anti-IKZF4 rabbit polyclonal antibody

别名: EOS; ZNFN1A4

抗原: IKZF4

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Members of the Ikaros (ZNFN1A1; MIM 603023) family of transcription factors, which includes Eos, are expressed in lymphocytes and are implicated in the control of lymphoid development.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	IKZF4
<b>Immunogen:</b>	Synthetic peptide of human IKZF4
<b>Full name:</b>	IKAROS family zinc finger 4 (Eos)
<b>Synonyms :</b>	EOS; ZNFN1A4
<b>SwissProt:</b>	Q9H2S9
<b>ELISA Recommended diluti</b>	2000-5000

<b>on:</b>	
<b>IHC positive control:</b>	Human liver cancer and Human gastric cancer
<b>IHC Recommend dilution:</b>	20-100

